

Anti-Anthrax Protective gen antibody [5C1-C8-H8] (STJA0032497)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	B.anthraxis

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	ELISA 1:10000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anthrax Protective Antigen antibody standard curve mouse STJA0032497