

## Anti-Lamin B2 antibody [8F6-E8-F12] (STJA0032503)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Human/Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">84823</a>
Gene Symbol	<a href="#">LMNB2</a>
UniProt ID	<a href="#">LMNB2_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Lamin B2 antibody western blot HCT116 mouse embryo Jurkat STJA0032503



Lamin B2 antibody immunoprecipitation HeLa lysates STJA0032503