

## Anti-Hsp90 alpha antibody [6H7-A3-2G6] (STJA0032512)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Monkey

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.6 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3320</a>
Gene Symbol	<a href="#">HSP90AA1</a>
UniProt ID	<a href="#">HS90A_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Hsp90 alpha antibody western blot rat brain Jurkat  
CHO-K1 COS7 C6 3T3 HeLa STJA0032512