

Anti-Ado-trastuzumab emtansine-HER2 antibody (STJA0032514)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FUNC
Host / Source	Humanized
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	0.01M PBS, pH 7.2
Isotype	IgG1
Storage Instruction	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at -80°C. Avoid repeated freeze/thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---