

## Anti-Lgi4 antibody [11H12] (STJA0032585)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ICC/IF          |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 1 mg/mL                                                                                  |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IF 1:50-1:200                                                         |
| Formulation         | Liquid in PBS containing 0.09% sodium azide as a preservative.                           |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">243914</a>     |
| Gene Symbol | <a href="#">Lgi4</a>       |
| UniProt ID  | <a href="#">LGI4_MOUSE</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Lgi4 antibody western blot COS7 lysates  
 STJA0032585