

**Anti-PUS7 antibody [4C6] (STJA0032686)****GENERAL INFORMATION**

|               |                        |
|---------------|------------------------|
| Product Type  | Primary antibodies     |
| Applications  | WB/FC                  |
| Host / Source | Mouse                  |
| Reactivity    | Human/Monkey/Mouse/Rat |

**PRODUCT PROPERTIES**

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.6 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>FC 1:50-1:100                                                         |
| Formulation         | Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.            |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">54517</a>      |
| Gene Symbol | <a href="#">PUS7</a>       |
| UniProt ID  | <a href="#">PUS7_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|