

Anti-VSV G tag-Tag antibody [4A10-9A6-6H5] (STJA0032754)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IP/ICC/IF
Host / Source	Mouse
Reactivity	Species-independent

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000 IF 1:100-1:500 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



VSV G tag antibody western blot analysis
STJA0032754



VSV G tag antibody immunofluorescence 293T
cells 4A10 STJA0032754



VSV-G antibody immunoprecipitation 293 cells
STJA0032754