

Anti-CD23 antibody [3D1-2G4-9B6] (STJA0032768)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IF/IHC-F/IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	2208
Gene Symbol	FCER2
UniProt ID	FCER2_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD23 antibody immunofluorescence human stomach STJA0032768




CD23 antibody immunofluorescence rat lung tissue 1E9 STJA0032768




CD23 antibody immunohistochemistry human lung tissue STJA0032768



CD23 antibody immunohistochemistry rat kidney tissue STJA0032768



CD23 antibody immunohistochemistry mouse lung tissue STJA0032768



CD23 antibody immunohistochemistry human tonsil tissue STJA0032768