

## Anti-CD21 antibody [3B1-5G10-6H1] (STJA0032776)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P/ICC/IF       |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 1 mg/mL                                                                                  |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | IHC 1:50-1:100<br>IF 1:50-1:200                                                          |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                           |
|-------------|---------------------------|
| Gene ID     | <a href="#">1380</a>      |
| Gene Symbol | <a href="#">CR2</a>       |
| UniProt ID  | <a href="#">CR2_HUMAN</a> |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD21 antibody immunohistochemistry human stomach tissue STJA0032776




CD21 antibody immunohistochemistry human tonsils STJA0032776



CD21 antibody immunohistochemistry mouse kidney tissue STJA0032776



CD21 antibody immunofluorescence mouse liver tissue 2C5 STJA0032776



CD21 antibody immunohistochemistry human tonsil tissue STJA0032776