

Anti-NSE antibody [4C7-4F10-10F9] (STJA0032779)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	2026
Gene Symbol	ENO2
UniProt ID	ENOG_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Enolase2 antibody western blot COS7 C6 HeLa STJA0032779



NSE antibody immunohistochemistry human lung small cell carcinoma tissue STJA0032779



Enolase 4C7 antibody western blot CHO-K1 rat brain C6 3T3 HeLa STJA0032779



NSE 4C7 antibody immunocytochemistry HeLa cells STJA0032779