

Anti-MAP2 antibody [3B5-9D6-10E6] (STJA0032780)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	4133
Gene Symbol	MAP2
UniProt ID	MTAP2_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




MAP2 antibody immunohistochemistry human tonsil tissue STJA0032780




MAP2 antibody immunohistochemistry human tonsils STJA0032780




MAP2 antibody immunohistochemistry mouse brain tissue STJA0032780



MAP2 antibody immunofluorescence mouse brain tissue 7D4 STJA0032780



MAP2 antibody immunofluorescence rat brain STJA0032780



MAP2 antibody immunohistochemistry human brain tissue STJA0032780