

Anti-Cytokeratin 17 antibody [8C6-4H2-5C1] (STJA0032784)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	3872
Gene Symbol	KRT17
UniProt ID	K1C17_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Cytokeratin 17 antibody immunohistochemistry human liver tissue STJA0032784



CK17 antibody immunofluorescence human breast tissue STJA0032784



Cytokeratin 17 antibody western blot HeLa MCF-7 293T STJA0032784



Cytokeratin 17 antibody immunoprecipitation HeLa lysates STJA0032784