

Anti-CD15 antibody [8C10-6F3-2F5] (STJA0032792)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	2526
Gene Symbol	FUT4
UniProt ID	FUT4_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD15 antibody immunohistochemistry human lung cancer tissue STJA0032792




CD15 antibody immunofluorescence human liver cancer tissue STJA0032792




CD15 antibody immunohistochemistry rat brain tissue STJA0032792



CD15 antibody immunohistochemistry human stomach 8C10 STJA0032792



CD15 antibody immunohistochemistry human stomach 8C10 STJA0032792



CD15 antibody immunohistochemistry human stomach 8C10 STJA0032792