

Anti-Hydroxyacid Oxidase 1 antibody [4B4-3F7-9A6] (STJA0032829)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	54363
Gene Symbol	HAO1
UniProt ID	HAOX1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Hydroxyacid Oxidase 1 antibody immunohistochemistry human liver STJA0032829




Hydroxyacid Oxidase 1 antibody immunohistochemistry human tonsils STJA0032829




HAO1 antibody immunohistochemistry mouse heart tissue STJA0032829



Hydroxyacid Oxidase 1 antibody immunofluorescence human appendix STJA0032829



HAO1 antibody immunofluorescence mouse spleen STJA0032829



HAO1 antibody western blot mouse rat liver STJA0032829