

Anti-Cytochrome C antibody [7C10-3H3-9A3] (STJA0032831)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Chicken

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	54205
Gene Symbol	CYCS
UniProt ID	CYC_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Cytochrome c antibody western blot T47D 3T3 COS7 C6 HeLa STJA0032831



Cytochrome C antibody western blot HeLa T47D K562 mouse rat brain STJA0032831



Cytochrome C 7C10-3H3-9A3 antibody immunocytochemistry HeLa STJA0032831



Cytochrome C antibody immunohistochemistry human tonsil tissue STJA0032831