

Anti-Flag-Tag antibody [7E1-10D2-3H1] (STJA0032833)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/ICC/IF/IP
Host / Source	Mouse
Reactivity	Species-independent

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000 IF 1:100-1:500 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Flag Tag 7E1 antibody immunofluorescence 293 cells STJA0032833



Flag Tag antibody immunoprecipitation 293 cells STJA0032833



Flag Tag antibody western blot CHO-K1 cells STJA0032833