

Anti-beta Actin antibody [10F1-3A3-3B1] (STJA0032843)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Rat/Mouse/Monkey/Dog/Chicken/Hamster/Rabbit/Insect

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.76 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:5000-1:10000 IHC 1:50-1:200 IF 1:100-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	60
Gene Symbol	ACTB
UniProt ID	ACTB_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




beta actin antibody western blot MCF-7 3T3 rat brain Jurkat HeLa STJA0032843



beta Actin antibody immunohistochemistry human ovary STJA0032843



Beta-Actin antibody western blot HeLa rat mouse chicken rabbit STJA0032843



BetaActin antibody western blot CHO-K1 rat brain C6 3T3 HeLa STJA0032843