

## Anti-GAPDH antibody [8A3-6C3-8D6] (STJA0032844)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/ICC/IF/IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse/Monkey/Dog/Chicken/Hamster/Rabbit/Pig/Sheep/Insect/Yeast

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:5000-1:10000 IHC 1:50-1:200 IF 1:100-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2597</a>
Gene Symbol	<a href="#">GAPDH</a>
UniProt ID	<a href="#">G3P_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

