

Anti-beta I Tubulin antibody [1H10-8B8-3D2] (STJA0032848)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	81027
Gene Symbol	TUBB1
UniProt ID	TBB1_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




beta I Tubulin antibody immunohistochemistry human liver cancer STJA0032848




beta I Tubulin antibody western blot CHO-K1 Jurkat 3T3 C6 HeLa STJA0032848



beta I Tubulin antibody western blot HeLa mouse rat brain STJA0032848



beta I Tubulin antibody western blot A549 K562 THP-1 STJA0032848



beta I Tubulin antibody immunocytochemistry HeLa cells STJA0032848