

## Anti-alpha Tubulin antibody [4C7-3B1-10G4] (STJA0032860)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Mouse              |
| Reactivity    | Zebrafish          |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 1 mg/mL                                                                                  |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000                                                                          |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">10376</a>       |
| Gene Symbol | <a href="#">TUBA1B</a>      |
| UniProt ID  | <a href="#">TBA1B_HUMAN</a> |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



alpha tubulin antibody western blot zebrafish skeletal muscle STJA0032860



alpha tubulin antibody western blot Raji MCF-7 293T HCT116 HeLa STJA0032860