

Anti-Phospho-Phospho-Tyrosine antibody [6G6-3D7-7A7] (STJA0032875)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Mouse
Reactivity	Transfected

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-tyrosine antibody western blot Jurkat 293T rat mouse brain STJA0032875



Phospho-tyrosine antibody immunohistochemistry human breast carcinoma STJA0032875