

Anti-CD1a antibody [4C3-9C8-6A10] (STJA0032878)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	909
Gene Symbol	CD1A
UniProt ID	CD1A_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD1a antibody immunohistochemistry human tonsil tissue STJA0032878




CD1a antibody immunohistochemistry human tonsils STJA0032878



CD1a antibody immunohistochemistry mouse heart tissue STJA0032878



CD1a antibody immunofluorescence mouse heart tissue 9H6 STJA0032878



CD1a antibody immunohistochemistry human tonsils STJA0032878