

Anti-CD4 antibody [7H9-1H6-7C6] (STJA0032880)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	920
Gene Symbol	CD4
UniProt ID	CD4_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD4 antibody immunohistochemistry human stomach tissue STJA0032880




CD4 antibody immunohistochemistry human tonsils STJA0032880




CD4 antibody immunohistochemistry mouse brain tissue STJA0032880



CD4 antibody immunofluorescence mouse colon tissue 11A1 STJA0032880



CD4 antibody immunofluorescence rat lung 11A1 STJA0032880



CD4 antibody immunohistochemistry human amygdala STJA0032880