

Anti-CD5 antibody [10F9-5E2-3D3] (STJA0032881)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	921
Gene Symbol	CD5
UniProt ID	CD5_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD5 antibody immunohistochemistry human liver cancer tissue STJA0032881




CD5 antibody immunohistochemistry human tonsils STJA0032881



CD5 antibody immunohistochemistry mouse kidney STJA0032881



CD5 antibody immunofluorescence mouse heart tissue 10G8 STJA0032881



CD5 antibody immunohistochemistry human tonsil carcinoma STJA0032881