

Anti-Phospho-CaMKII-Thr287 antibody [4A8-3A1-10C8] (STJA0032887)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID
[818](#)
[816](#)
[817](#)
[CAMK2G](#)
[CAMK2B](#)
[KCC2G_HUMAN](#)
[KCC2B_HUMAN](#)
[KCC2D_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-CaMKII antibody immunohistochemistry
 rat brain STJA0032887



Phospho-CaMKII antibody immunohistochemistry
 human tonsils STJA0032887