

## Anti-Phospho-ERK1/2-Tyr222/Tyr205 antibody [1H4-6D7-9D8] (STJA0032903)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">5594</a> <a href="#">5595</a>
Gene Symbol	<a href="#">MAPK1</a> <a href="#">MAPK3</a>
UniProt ID	<a href="#">MK01_HUMAN</a> <a href="#">MK03_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-ERK1/2 antibody immunohistochemistry human colon carcinoma STJA0032903



Phospho-ERK1/2 antibody immunohistochemistry human tonsils STJA0032903