

## Anti-GSK3 beta antibody [9F9-5G6-7G2] (STJA0032904)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P           |
| Host / Source | Mouse              |
| Reactivity    | Human/Rat/Mouse    |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 1 mg/mL                                                                                  |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100                                                        |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">2932</a>        |
| Gene Symbol | <a href="#">GSK3B</a>       |
| UniProt ID  | <a href="#">GSK3B_HUMAN</a> |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




GSK3 beta antibody immunohistochemistry human breast carcinoma STJA0032904



GSK3 beta antibody immunohistochemistry human tonsils STJA0032904



GSK3 beta antibody western blot goat muscle CHO-K1 C6 HeLa 3T3 STJA0032904



GSK3 beta antibody western blot HeLa A549 HL-60 U2OS STJA0032904