

Anti-DDIT3 antibody [7G7-3B2-2F10] (STJA0032921)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1649
Gene Symbol	DDIT3
UniProt ID	DDIT3_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




DDIT3 antibody immunohistochemistry human stomach tissue STJA0032921




DDIT3 antibody immunofluorescence mouse brain tissue STJA0032921



CHOP antibody immunohistochemistry human stomach carcinoma tissue STJA0032921



CHOP antibody immunohistochemistry human pancreatic carcinoma STJA0032921



DDIT3 antibody western blot mouse liver STJA0032921