

Anti-HP1 alpha antibody [1E8-5C3-10E1] (STJA0032938)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	23468
Gene Symbol	CBX5
UniProt ID	CBX5_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




HP1 alpha antibody immunohistochemistry human colon carcinoma STJA0032938




HP1 alpha antibody immunohistochemistry human tonsils STJA0032938



HP1 alpha antibody immunohistochemistry human placenta STJA0032938



HP1 alpha antibody western blot COS7 HeLa STJA0032938



HP1-alpha antibody western blot HeLa A549 HL-60 U2OS C6 STJA0032938