

Anti-Hsp90 alpha antibody [5D6-2G8-6D5] (STJA0032960)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	3320
Gene Symbol	HSP90AA1
UniProt ID	HS90A_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Hsp90 alpha antibody immunohistochemistry human colon carcinoma STJA0032960



Hsp90 alpha antibody western blot K562 Ramos Jurkat rat brain STJA0032960



Hsp90 alpha antibody western blot HeLa A549 HL60 U2OS C6 STJA0032960