

Anti-PI3 Kinase p85 alpha antibody [1C8-2H4-9D9] (STJA0032965)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5295
Gene Symbol	PIK3R1
UniProt ID	P85A_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




PI3 Kinase p85 alpha antibody immunohistochemistry rat liver STJA0032965



PI3 Kinase p85 alpha antibody immunohistochemistry human tonsils STJA0032965



PI3 Kinase p85 alpha antibody western blot Jurkat rat brain Raw2647 STJA0032965



PI3 Kinase p85 alpha antibody western blot HeLa A549 HL-60 U2OS C6 STJA0032965