

Anti-Galectin 3 antibody [7G9-1C9-7B4] (STJA0032978)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 ELISA 1:10000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	3958
Gene Symbol	LGALS3
UniProt ID	LEG3_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Galectin 3 antibody immunohistochemistry human colon STJA0032978



Galectin 3 antibody immunohistochemistry human tonsils STJA0032978



Galectin 3 antibody western blot CHO-K1 HeLa STJA0032978



Galectin 3 antibody western blot MCF-7 3T3 STJA0032978