

Anti-Ki67 antibody [12B2] (STJA0032986)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/ICC/IF
Host / Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IHC 1:50-1:100 IF 1:50-1:200
Formulation	Rat IgG in phosphate buffered saline, pH 7.4, 150mM NaCl and 0.02% sodium azide.
Isotype	IgG2b
Storage Instruction	Store at 2-8°C. Do not freeze.

TARGET INFORMATION

Gene ID	4288
Gene Symbol	MKI67
UniProt ID	KI67_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Ki67 antibody immunohistochemistry human liver cancer Chimera mouse STJA0032986