

Anti-F4/80 antibody [31B1] (STJA0032987)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-P/ICC/IF/FC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity Chromatography |
| Dilution Range | IHC 1:50-1:100 IF 1:50-1:200 FC 1:50-1:100 |
| Formulation | Rat IgG in phosphate buffered saline, pH 7.4, 150mM NaCl and 0.02% sodium azide. |
| Isotype | IgG2a |
| Storage Instruction | Store at 2-8°C. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 13733 |
| Gene Symbol | Adgre1 |
| UniProt ID | AGRE1_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|