

Anti-BrdU antibody [41B2] (STJA0032988)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IF
Host / Source	Mouse
Reactivity	All

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IF 1:50-1:200
Formulation	Mouse IgG in phosphate buffered saline, pH 7.4, 150mM NaCl and 0.02% sodium azide.
Isotype	IgG1
Storage Instruction	Store at 2-8°C. Do not freeze.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---