

## Anti-Flag-Tag antibody [10E1-1D8-5F5] (STJA0032989)

### GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Species-independent

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000 IP 1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Flag tag antibody western blot CHO-K1 transfected STJA0032989



Flag tag antibody western blot CHO-K1 cells STJA0032989