

Anti-CD44 antibody [Hermes-3] (STJA0033002)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/WB/ICC/IF/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 FC 1:50-1:100
Formulation	Liquid in PBS containing 0.09% sodium azide, pH 7.3.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	960
Gene Symbol	CD44
UniProt ID	CD44_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD44 antibody western blot A549 HeLa
STJA0033002



CD44 antibody immunohistochemistry human lung adenocarcinoma tissue
STJA0033002



CD44 antibody immunohistochemistry human skin tissue
STJA0033002



CD44 antibody immunofluorescence HeLa cells
STJA0033002