

Anti-Rac1 antibody [1F4-5C8-6D4] (STJA0033004)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 FC 1:50-1:100
Formulation	Liquid in PBS containing 0.09% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5879
Gene Symbol	RAC1
UniProt ID	RAC1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




RAC1 antibody western blot A431 NIH/3T3 C6 STJA0033004



RAC1 antibody immunohistochemistry human skin STJA0033004



RAC1 antibody flow cytometry U87 MG cells STJA0033004



RAC1 antibody immunohistochemistry human stomach STJA0033004