

Anti-MonoMethyl-Histone H3-Lys27 antibody [4K2-P7-G3] (STJA0033024)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	All

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000
Formulation	Mouse IgG in phosphate buffered saline, glycerol and BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---