

Anti-MonoMethyl-Histone H2B-Lys15 antibody [2Y1-Z2-D4] (STJA0033044)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000
Formulation	Mouse IgG in phosphate buffered saline, glycerol and BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	3017/8339/8343/8344/8346/8347
Gene Symbol	H2BC4.H2BC6.H2BC7.H2BC8.H2BC10
UniProt ID	H2B1C_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---