

Anti-IL-1 beta antibody [R09-4F-7] (STJA0033050)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG1k
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	3553
Gene Symbol	IL1B
UniProt ID	IL1B_HUMAN

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



IL-1 beta antibody western blot THP-1 cells
 STJA0033050