

Anti-NeuN antibody [R07-8O-7] (STJA0033053)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC-P 1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG1k
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	146713
Gene Symbol	RBFOX3
UniProt ID	RFOX3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



NeuN antibody immunohistochemistry mouse brain STJA0033053