

Anti-BRG1 antibody [R05-9P-8] (STJA0033059)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2Bk
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	6597
Gene Symbol	SMARCA4
UniProt ID	SMCA4_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



BRG1 antibody western blot HeLa HepG2 C6
 Jurkat STJA0033059