

Anti-GATA3 antibody [R01-7I-9] (STJA0033063)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC-P 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2bk
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	2625
Gene Symbol	GATA3
UniProt ID	GATA3_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



GATA3 antibody immunohistochemistry mouse kidney STJA0033063