

Anti-CCR4 antibody [R07-5I-7] (STJA0033070)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	FC 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2bk
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	1233
Gene Symbol	CCR4
UniProt ID	CCR4_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CCR4 antibody flow cytometry PBMC negative control STJA0033070