

## Anti-ILK antibody [R02-9F-3] (STJA0033073)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2Ak
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3611</a>
Gene Symbol	<a href="#">ILK</a>
UniProt ID	<a href="#">ILK_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



ILK antibody western blot HeLa SH-SY5Y 4T1 RAW2647 STJA0033073