

## Anti-BACH1 antibody [R05-1W-6] (STJA0033079)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG1k
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">571</a>
Gene Symbol	<a href="#">BACH1</a>
UniProt ID	<a href="#">BACH1_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



BACH1 antibody western blot HeLa Jurkat HEK-293T STJA0033079