

Anti-AKR7A2 antibody [R03-5C-3] (STJA0033082)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ICC/IF
	Host / Source	Mouse
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein G
	Dilution Range	WB 1:500-1:1000 ICC:IF 1:50
	Formulation	Liquid in PBS, Glycerol and BSA
	Isotype	IgG1k
	Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.
TARGET INFORMATION		
	Gene ID	8574
	Gene Symbol	AKR7A2
	UniProt ID	ARK72_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.