

Anti-Glycophorin A antibody [R07-2E-1] (STJA0033086)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC-P 1:50-1:100 FC 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG1k
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	2993
Gene Symbol	GYPA
UniProt ID	GLPA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---