

Anti-CEP55 antibody [R07-4J-2] (STJA0033124)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB/ICC/IF |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein G |
| Dilution Range | WB 1:500-1:1000 ICC:IF 1:50 |
| Formulation | Liquid in PBS, Glycerol and BSA |
| Isotype | IgG1k |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 55165 |
| Gene Symbol | CEP55 |
| UniProt ID | CEP55_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|